

## Product datasheet for **SM3151P**

### IL2 Receptor alpha (IL2RA) Mouse Monoclonal Antibody [Clone ID: MEM-140]

#### Product data:

Product Type:	Primary Antibodies
Clone Name:	MEM-140
Applications:	FC, IP
Recommended Dilution:	<b>Immunoprecipitation:</b> excellent for immunoisolation of CD25+cells <b>Flow Cytometry:</b> 1 µg/ml.
Reactivity:	Human
Host:	Mouse
Isotype:	IgM
Clonality:	Monoclonal
Immunogen:	PHA-activated peripheral blood leucocytes
Specificity:	The antibody MEM-140 reacts with CD25 (Interleukin-2 receptor a chain), a 55 kDa type I transmembrane glycoprotein expressed on activated B and T lymphocytes, activated monocytes/macrophages and on CD4+ T lymphocytes (T regulatory cells); it is lost on resting B and T lymphocytes.
Formulation:	Tris buffered saline (TBS) with 15 mM sodium azide, approx. pH 8.0 State: Purified State: Liquid Ig fraction
Concentration:	lot specific
Purification:	Thiophilic adsorption-affinity chromatography; purity: > 95% (by SDS-PAGE)
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>
Stability:	Shelf life: one year from despatch.
Gene Name:	interleukin 2 receptor subunit alpha
Database Link:	<a href="#">Entrez Gene 3559 Human P01589</a>




[View online »](#)

**Background:**

CD25 (IL2Ralpha, Tac) is a ligand-binding alpha subunit of interleukin 2 receptor (IL2R). Together with beta and gamma subunit CD25 constitutes the high affinity IL2R, whereas CD25 alone serves as the low affinity IL2R. CD25 expression rapidly increases upon T cell activation. The 55 kDa CD25 molecule is enzymatically cleaved and shed from the cell surface as a soluble 45 kDa s-Tac, whose concentration in serum can be used as a marker of T cell activation. Expression of CD25 indicates the neoplastic phenotype of mast cells. Humanized anti CD25 antibodies represent a useful tool to reduce the incidence of allograft rejection as well as the severity of graft versus host reaction, and radioimmunoconjugates of anti-CD25 antibodies can be used against CD25 expressing lymphomas.

**Synonyms:**

Interleukin-2 receptor alpha chain, IL-2 receptor alpha subunit, IL-2-RA, IL2-RA, p55, TAC antigen

**Product images:**

Surface staining of CD25 in PHA activated PBMC with anti-CD25 (MEM-140) purified, GAM-APC.