

Product datasheet for **SM3157P**

Human IgM (Mu chain, Fab-Region) Mouse Monoclonal Antibody [Clone ID: MA2]

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgM (Mu chain, Fab-Region) Mouse Monoclonal Antibody [Clone ID: MA2]
Clone Name:	MA2
Applications:	ELISA
Recommended Dilution:	ELISA.
Reactivity:	Human
Host:	Mouse
Immunogen:	IgM from Human serum.
Isotype:	IgG2b
Formulation:	Tris buffered saline (TBS) with 15 mM Sodium Azide, approx. pH 8.0 State: Purified State: Liquid purified Ig fraction (> 95% pure by SDS-PAGE).
Concentration:	lot specific
Purification:	Precipitation methods and Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)