

Product datasheet for **SM427**

MHC Class II (monomorphic) Mouse Monoclonal Antibody [Clone ID: vpg3]

Product data:

Product Type:	Primary Antibodies
Clone Name:	vpg3
Applications:	FC, IHC, IP
Recommended Dilution:	Flow Cytometry. Immunoprecipitation. Immunohistochemistry on frozen sections.
Reactivity:	Feline
Host:	Mouse
Isotype:	IgG1
Clonality:	Monoclonal
Immunogen:	IL-2 dependent feline T cells. Spleen cells from immunised BALB/c mice were fused with cells of the NSO mouse myeloma cell line.
Specificity:	This antibody recognises MHC Class II antigens.
Formulation:	Containing 0.2M Tris/HCl pH7.4 and 5-10% foetal calf serum and 0.09% Sodium Azide State: Supernatant State: Liquid Tissue Culture Supernatant
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.
Synonyms:	MHC Class 2, HLA Class II



[View online »](#)