


Product datasheet for SO25000-100

Tyrosine Hydroxylase (TH) Sheep Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IF, IHC, WB
Reactivity:	Human, Mouse, Primate, Rat
Host:	Sheep
Isotype:	IgG
Clonality:	Polyclonal
Formulation:	State: Aff - Purified
Gene Name:	tyrosine hydroxylase
Database Link:	Entrez Gene 7054 Human
Synonyms:	Tyrosine 3-hydroxylase, TYH


Product images:




Detection of tyrosine hydroxylase in PC12 cell lysate by Western Blot. Band is at approximately 60 kDa.




[View online »](#)



TH staining of neurons in the mouse SNc.



Detection of tyrosine hydroxylase in PC12 cell lysate by Western Blot. Band is at approximately 60 kDa.



TH staining of neurons in the mouse SNc.