

## Product datasheet for **SP1018HRP**

### Mouse IgM (Mu chain specific), adsorbed Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgM (Mu chain specific), adsorbed Goat Polyclonal Antibody
Applications:	ELISA, IHC
Recommended Dilution:	<b>ELISA</b> (1/1000). <b>Immunohistochemistry on Frozen and Paraffin Tissues</b> (1/50).
Reactivity:	Mouse
Host:	Goat
Immunogen:	Monoclonal Mouse IgM
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction Stabilizer: HRP stabilizer Preservative: 0.01% Thimerosal Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C. <b>DO NOT FREEZE!</b>



[View online »](#)