

Product datasheet for **SP1020BS**

Rat IgG (H+L chain, adsorbed) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG (H+L chain, adsorbed) Goat Polyclonal Antibody
Applications:	IF, IHC, WB
Recommended Dilution:	ELISA: 1/5000-1/20000. Western Blot. Immunofluorescence. Immunohistochemistry on Frozen Sections. Immunohistochemistry on Paraffin Sections.
Reactivity:	Rat
Host:	Goat
Formulation:	PBS Label: Biotin State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Note:	Antiserum Preparation: Antisera to rat IgG were raised by repeated immunisation of goats with purified antigen.



[View online »](#)