

## Product datasheet for **SP1020HRP**

### Rat IgG (H+L chain, adsorbed) Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG (H+L chain, adsorbed) Goat Polyclonal Antibody
Applications:	ELISA, IHC
Recommended Dilution:	<b>ELISA:</b> 1/1000. <b>Immunohistochemistry on Frozen Sections:</b> 1/50. <b>Immunohistochemistry on Paraffin Sections.</b>
Reactivity:	Rat
Host:	Goat
Immunogen:	Rat IgG
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction Stabilizer: 1% BSA Preservative: 0.01% Thiomersal, HRP Stabilizer Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Note:	<b>Antiserum Preparation:</b> Antisera to rat IgG were raised by repeated immunisation of goats with highly purified antigen.



[View online »](#)