

Product datasheet for **SP1021HRP**

Rat IgG (H+L chain), Fab Fragment Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG (H+L chain), Fab Fragment Rabbit Polyclonal Antibody
Recommended Dilution:	ELISA: 1/1000 - 1/2000. Immunohistochemistry on frozen sections: 1/50. Immunohistochemistry on paraffin sections: 1/50. Immunohistochemistry on resin sections: 1/50.
Reactivity:	Rat
Host:	Rabbit
Immunogen:	Purified rat immunoglobulin G
Formulation:	PBS, pH 7.2 containing 0.01% Thiomersal and a HRP stabiliser Label: HRP State: Liquid Fab fragment of IgG Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)