

Product datasheet for **SP1021P**

Rat IgG (H+L chain), Fab Fragment Rabbit Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Rat IgG (H+L chain), Fab Fragment Rabbit Polyclonal Antibody
Recommended Dilution:	ELISA: 1/250 - 1/1000. Immunodiffusion.
Reactivity:	Rat
Host:	Rabbit
Immunogen:	Purified polyclonal rat IgG
Formulation:	PBS, pH 7.4 containing 0.09% Sodium Azide State: Purified State: Liquid Fab fragment of IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)