

# Product datasheet for SP1025B

## Hamster IgG (adsorbed) Goat Polyclonal Antibody

### **Product data:**

Product Type:	Secondary Antibodies
Product Name:	Hamster IgG (adsorbed) Goat Polyclonal Antibody
Applications:	ELISA, FC
<b>Recommended Dilution:</b>	ELISA: 1/1000-1/5000.
	<b>Flow Cytometry:</b> Use 50 μl of 1/25-1/50 diluted antibody to label 1 <b>6</b> cells in 100 μl.
Reactivity:	Hamster
Host:	Goat
Immunogen:	Hamster IgG
lsotype:	lgG
Formulation:	PBS Label: Biotin State: Liquid purified IgG fraction Stabilizer: 1% BSA Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.


#### OriGene Technologies, Inc.


9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US

### **Product images:**





This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2022 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US