

## Product datasheet for **SP1025F**

### Hamster IgG (adsorbed) Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Hamster IgG (adsorbed) Goat Polyclonal Antibody
Applications:	FC
Recommended Dilution:	Flow Cytometry: 1/50; Use 50µl of the suggested working dilution to label 1 x 10e6 cells in 100µl.
Reactivity:	Hamster
Host:	Goat
Immunogen:	Hamster IgG
Formulation:	PBS containing 0.09% Sodium Azide and 1% Bovine Serum Albumin Label: FITC State: Liquid purified IgG Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)