

Product datasheet for **SP1025HRP**

Hamster IgG (adsorbed) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Hamster IgG (adsorbed) Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/1000 - 1/5000.
Reactivity:	Hamster
Host:	Goat
Immunogen:	Highly purified Hamster IgG
Formulation:	PBS containing 1% BSA as stabilizer and 0.01% thiomersal as preservative Label: HRP State: Liquid purified IgG fraction Label: Horseradish peroxidase
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)