

Product datasheet for **SP1032P**

Human IgG (Fc specific) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Human IgG (Fc specific) Goat Polyclonal Antibody
Applications:	ELISA, IHC
Recommended Dilution:	ELISA: 1-5 µg/ml. Immunohistochemistry on Frozen Sections.
Reactivity:	Human
Host:	Goat
Immunogen:	Human IgG
Formulation:	PBS, pH 7.4 State: Aff - Purified State: Liquid purified IgG fraction from whole serum Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)