

Product datasheet for **SP1041P**

Mouse IgG (Rat adsorbed) Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Mouse IgG (Rat adsorbed) Goat Polyclonal Antibody
Applications:	ELISA, FC
Recommended Dilution:	ELISA. Flow Cytometry.
Reactivity:	Mouse
Host:	Goat
Immunogen:	Mouse IgG
Formulation:	PBS containing 0.09% Sodium Azide State: Aff - Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatograph
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)