

Product datasheet for **SP1047APX**

Bovine IgG Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgG Sheep Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/5000-1/50000. Western Blot: 1/5000-1/50000.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified Bovine IgG.
Isotype:	IgG
Formulation:	50mM HEPES, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂
	Label: AP State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide Label: Alkaline Phosphatase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	AP
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!
Note:	Antiserum Preparation: Antisera to bovine IgG were raised by repeated immunisation of sheep with highly purified antigen.



[View online »](#)