

Product datasheet for **SP1047F**

Bovine IgG Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgG Sheep Polyclonal Antibody
Applications:	FC, IHC
Recommended Dilution:	Flow Cytometry. Immunohistochemistry on Frozen Sections: 1/200-1/2,000.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified Bovine IgG.
Isotype:	IgG
Formulation:	PBS Label: FITC State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)