

## Product datasheet for **SP1047HRP**

### Bovine IgG Sheep Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgG Sheep Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	<b>ELISA:</b> 1/10,000-1/100,000. <b>Western Blot.</b>
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified Bovine IgG.
Isotype:	IgG
Formulation:	PBS containing 0.1% ProClin™300 Label: HRP State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>



[View online »](#)