

Product datasheet for **SP1047P**

Bovine IgG Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgG Sheep Polyclonal Antibody
Applications:	ELISA, ID
Recommended Dilution:	ELISA: 1/100-1/10,000. Immunodiffusion.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified Bovine IgG.
Isotype:	IgG
Formulation:	PBS State: Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)