

Product datasheet for **SP1048APX**

Bovine IgM Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgM Sheep Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/1000-1/10000. Western Blot: 1/5000-1/50000.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified Bovine IgM.
Isotype:	IgG
Formulation:	50mM HEPES, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂
	Label: AP State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide Label: Alkaline Phosphatase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	AP
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)