

## Product datasheet for **SP1048F**

### Bovine IgM Sheep Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgM Sheep Polyclonal Antibody
Applications:	FC, IHC
Recommended Dilution:	<b>Flow Cytometry:</b> 1/20-1/100. <b>Immunohistochemistry on frozen sections:</b> 1/200-1/2000.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified Bovine IgM
Formulation:	PBS containing 0.09% Sodium Azide as preservative and 0.2% BSA as stabilizer Label: FITC State: Liquid IgG fraction Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)