

Product datasheet for **SP1048HRP**

Bovine IgM Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgM Sheep Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10,000 - 1/100,000.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified bovine IgM
Formulation:	PBS containing 0.1% Proclin 300 and 0.2% Bovine Serum Albumin Label: HRP State: Liquid purified IgG Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)