

Product datasheet for **SP1048P**

Bovine IgM Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgM Sheep Polyclonal Antibody
Applications:	ELISA, ID
Recommended Dilution:	ELISA: 1/100-1/1000. Immunodiffusion.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified Bovine IgM
Formulation:	PBS, pH 7.2 containing 0.09% Sodium Azide as preservative State: Aff - Purified State: Liquid purified Ig fraction
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)