

Product datasheet for **SP1049HRP**

Bovine IgG1 Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgG1 Sheep Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10000-1/100000.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified Bovine IgG1
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction Preservative: 0.1% ProClin300 Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)