

Product datasheet for **SP1220**

Leumorphin Rabbit Polyclonal Antibody

Product data:

Product Type:	Primary Antibodies
Applications:	IHC
Recommended Dilution:	Immunohistochemistry on frozen sections: 1/200 - 1/2000.
Reactivity:	Human, Rat
Host:	Rabbit
Clonality:	Polyclonal
Immunogen:	Leumorphin
Specificity:	This antibody blocks specifically with leumorphin (rat only). Cells are found in the paraventricular and supraoptic nucleus and fibres in the lateral hypothalamus in rat brain of colchicine treated 4% FA perfused tissue.
Formulation:	PBS containing 0.09% Sodium Azide and 1% Bovine Serum Albumin State: Serum State: Liquid serum
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.
Stability:	Shelf life: one year from despatch.



[View online »](#)