

Product datasheet for **SP2000APX**

Bovine IgG2 Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgG2 Sheep Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/1000-1/10000.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified Bovine IgG2.
Isotype:	IgG
Formulation:	50mM HEPES, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂ Label: AP State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide Label: Alkaline Phosphatase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	AP
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)