

## Product datasheet for **SP2000F**

### Bovine IgG2 Sheep Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgG2 Sheep Polyclonal Antibody
Applications:	FC, IHC
Recommended Dilution:	Flow Cytometry: 1/20 - 1/100. Immunohistochemistry on frozen sections: 1/20 - 1/100.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified bovine IgG2
Formulation:	PBS containing 0.09% Sodium Azide Label: FITC State: Liquid Ig fraction Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	FITC
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing. This product is photosensitive and should be protected from light.



[View online »](#)