

Product datasheet for **SP2000HRP**

Bovine IgG2 Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgG2 Sheep Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1:10,000-1:100,000.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified bovine IgG2
Formulation:	Conjugated to Horseradish Peroxidase (HRP), in PBS with 0.1% ProClin 300 as stabiliser. Label: HRP State: Liquid
Concentration:	lot specific
Purification:	Purified IgG
Conjugation:	HRP
Storage:	Store the antibody at 4-8°C for one month or at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)