

Product datasheet for **SP2000P**

Bovine IgG2 Sheep Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Bovine IgG2 Sheep Polyclonal Antibody
Applications:	ELISA, ID
Recommended Dilution:	ELISA: 1/100 - 1/10000. Immunodiffusion.
Reactivity:	Bovine
Host:	Sheep
Immunogen:	Purified bovine IgG2
Formulation:	PBS containing 0.09% Sodium Azide State: Aff - Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)