

## Product datasheet for **SP2001HRP**

### Canine IgA Goat Polyclonal Antibody

#### Product data:

|                       |                                                                                                                                   |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Product Type:         | Secondary Antibodies                                                                                                              |
| Product Name:         | Canine IgA Goat Polyclonal Antibody                                                                                               |
| Applications:         | ELISA, IHC                                                                                                                        |
| Recommended Dilution: | ELISA: 1:10,000 - 1:80,000.<br>Immunohistochemistry on frozen sections: 1:200 - 1:500.                                            |
| Reactivity:           | Canine                                                                                                                            |
| Host:                 | Goat                                                                                                                              |
| Immunogen:            | Purified canine IgA                                                                                                               |
| Formulation:          | PBS, containing 0.1% ProClin 300 and 0.2% BSA<br>Label: HRP<br>State: Liquid purified IgG<br>Label: Horseradish Peroxidase        |
| Concentration:        | lot specific                                                                                                                      |
| Purification:         | Affinity chromatography                                                                                                           |
| Conjugation:          | HRP                                                                                                                               |
| Storage:              | Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.<br>Avoid repeated freezing and thawing. |



[View online »](#)