

Product datasheet for **SP2001P**

Canine IgA Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Canine IgA Goat Polyclonal Antibody
Applications:	ELISA, ID, IHC, WB
Recommended Dilution:	ELISA. Western Blot (Reducing conditions). Immunodiffusion. Immunohistochemistry on Frozen Sections: 1/200-1/500.
Reactivity:	Canine
Host:	Goat
Immunogen:	Purified Canine IgA.
Formulation:	PBS State: Aff - Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)