

## Product datasheet for **SP2004HRP**

### Canine IgM Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Canine IgM Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1/10,000-1/100,000.
Reactivity:	Canine
Host:	Goat
Immunogen:	Purified Canine IgM
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction. Stabilizer: 0.2% BSA Preservative: 0.1% Proclin™300 Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>



[View online »](#)