

Product datasheet for **SP2004P**

Canine IgM Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Canine IgM Goat Polyclonal Antibody
Applications:	ELISA, FC, ID, IHC
Recommended Dilution:	ELISA: 1/100-1/1000. Flow Cytometry. Immunodiffusion. Immunohistochemistry on Frozen Sections.
Reactivity:	Canine
Host:	Goat
Immunogen:	Purified Canine IgM.
Formulation:	PBS containing 0.09% Sodium Azide State: Aff - Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)