

Product datasheet for SP2006HRP

Equine IgA Goat Polyclonal Antibody

Product data:

Product Type: Secondary Antibodies

Product Name: Equine IgA Goat Polyclonal Antibody

Applications: ELISA

Recommended Dilution: ELISA: 1/10,000 - 1/100,000.

Reactivity: Equine Host: Goat

Immunogen: Purified equine IgA.

Formulation: PBS containing 0.1% ProclinTM300 and 0.2% Bovine Serum Albumin

Label: HRP

State: Liquid purified IgG Label: Horseradish Peroxidase

Concentration: lot specific

Purification: Affinity chromatography

Conjugation: HRP

Storage: Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer.

Avoid repeated freezing and thawing.



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com