

Product datasheet for **SP2006P**

Equine IgA Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Equine IgA Goat Polyclonal Antibody
Applications:	ELISA, ID
Recommended Dilution:	ELISA: 1/100 - 1/1000.
Reactivity:	Equine
Host:	Goat
Immunogen:	Purified equine IgA.
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)