

Product datasheet for **SP2007HRP**

Equine IgG Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Equine IgG Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10,000 - 1/100,000.
Reactivity:	Equine
Host:	Goat
Immunogen:	Purified equine IgG (T).
Formulation:	PBS containing 0.1% ProClin TM 300 and 0.2% Bovine Serum Albumin Label: HRP State: Liquid purified IgG Label: Horseradish Peroxidase
Reconstitution Method:	Restore with 1ml distilled water
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	HRP
Storage:	Prior to reconstitution store at 2-8°C. Following reconstitution store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)