

Product datasheet for **SP2012P**

Feline IgA Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Feline IgA Goat Polyclonal Antibody
Applications:	ELISA, ID
Recommended Dilution:	ELISA. Immunodiffusion.
Reactivity:	Feline
Host:	Goat
Immunogen:	Purified feline IgA.
Formulation:	PBS containing 0.09% Sodium Azide State: Aff - Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)