

## Product datasheet for **SP2013HRP**

### Feline IgG (Fc specific) Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Feline IgG (Fc specific) Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10,000.
Reactivity:	Feline
Host:	Goat
Immunogen:	Purified Feline IgG Fc fragment.
Formulation:	PBS containing 0.1% Proclin as preservative 300 and 0.2% BSA as stabilizer Label: HRP State: Liquid purified IgG fraction Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)