

Product datasheet for **SP2014F**

Feline IgM Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Feline IgM Goat Polyclonal Antibody
Applications:	FC, IHC
Recommended Dilution:	Flow Cytometry. Immunohistochemistry on Frozen Sections.
Reactivity:	Feline
Host:	Goat
Immunogen:	Purified Feline IgM
Formulation:	PBS containing 0.09% Sodium Azide as preservative and 0.2% BSA as stabilizer Label: FITC State: Liquid purified IgG fraction Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE! This product is photosensitive and should be protected from light.



[View online »](#)