

Product datasheet for **SP2014HRP**

Feline IgM Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Feline IgM Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10,000 - 1/100,000.
Reactivity:	Feline
Host:	Goat
Immunogen:	Purified Feline IgM
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.09% Proclin™300
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)