

Product datasheet for **SP2014P**

Feline IgM Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Feline IgM Goat Polyclonal Antibody
Applications:	ELISA, ID, IHC
Recommended Dilution:	ELISA.
Reactivity:	Feline
Host:	Goat
Immunogen:	Purified feline IgM
Formulation:	PBS containing 0.09% Sodium Azide State: Purified State: Liquid purified IgG
Concentration:	lot specific
Purification:	Affinity chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)