

Product datasheet for **SP2016APX**

Porcine IgA Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgA Goat Polyclonal Antibody
Applications:	ELISA, IHC, WB
Recommended Dilution:	ELISA: 1/1000-1/10000. Western Blot: 1/1000-1/10000. Immunohistochemistry on Frozen Sections: 1/100-1/1000.
Reactivity:	Porcine
Host:	Goat
Immunogen:	Purified Porcine IgA.
Isotype:	IgG
Formulation:	50mM HEPES, 0.1M NaCl, 1mM MgCl ₂ , 0.1mM ZnCl ₂ Label: AP State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide Label: Alkaline Phosphatase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	AP
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)