

Product datasheet for **SP2016B**

Porcine IgA Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgA Goat Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/10000-1/100000. Western Blot: 1/10000-1/100000
Reactivity:	Porcine
Host:	Goat
Immunogen:	Purified Porcine IgA
Isotype:	IgG
Formulation:	PBS Label: Biotin State: Liquid purified Ig fraction Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)