

## Product datasheet for **SP2016F**

### Porcine IgA Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgA Goat Polyclonal Antibody
Applications:	FC, IHC
Recommended Dilution:	<b>Flow Cytometry.</b> <b>Immunohistochemistry on Frozen Sections:</b> 1/200-1/2000.
Reactivity:	Porcine
Host:	Goat
Immunogen:	Purified Porcine IgA
Isotype:	IgG
Formulation:	PBS Label: FITC State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. This product is photosensitive and should be protected from light. Avoid repeated freezing and thawing.



[View online »](#)