

## Product datasheet for **SP2016HRP**

### Porcine IgA Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgA Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1/10,000-1/100,000.
Reactivity:	Porcine
Host:	Goat
Immunogen:	Purified Porcine IgA
Isotype:	IgG
Formulation:	PBS Label: HRP State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.1% Proclin 300 Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)