

## Product datasheet for **SP2016P**

### Porcine IgA Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgA Goat Polyclonal Antibody
Applications:	ELISA, ID
Recommended Dilution:	<b>ELISA:</b> 1/100-1/1000. <b>Immunodiffusion.</b>
Reactivity:	Porcine
Host:	Goat
Immunogen:	Purified Porcine IgA.
Isotype:	IgG
Formulation:	PBS State: Aff - Purified State: Liquid purified IgG fraction Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)