

## Product datasheet for **SP2017AP**

### Porcine IgG F(c) specific Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgG F(c) specific Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	<b>ELISA:</b> 1/5000-1/50000.
Reactivity:	Porcine
Host:	Goat
Immunogen:	Purified Porcine IgG
Formulation:	50mM HEPES, 0.1M NaCl, 1mM MgCl <sub>2</sub> , 0.1mM ZnCl <sub>2</sub> Label: AP State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide
Concentration:	Label: Alkaline Phosphatase lot specific
Purification:	Affinity Chromatography
Conjugation:	AP
Storage:	Store undiluted at 2-8°C. <b>DO NOT FREEZE!</b>



[View online »](#)