

Product datasheet for **SP2017B**

Porcine IgG F(c) specific Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgG F(c) specific Goat Polyclonal Antibody
Applications:	ELISA, WB
Recommended Dilution:	ELISA: 1/10000-1/100000. Western Blot: 1/10000-1/100000.
Reactivity:	Porcine
Host:	Goat
Immunogen:	Purified Porcine IgG.
Formulation:	PBS Label: Biotin State: Liquid purified IgG fraction Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Biotin
Storage:	Store undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)