

Product datasheet for **SP2017F**

Porcine IgG F(c) specific Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgG F(c) specific Goat Polyclonal Antibody
Applications:	FC, IHC
Recommended Dilution:	Flow Cytometry: 1/20-1/100. Immunohistochemistry on Frozen Sections: 1/20-1/100.
Reactivity:	Porcine
Host:	Goat
Immunogen:	Purified Porcine IgG.
Formulation:	PBS Label: FITC State: Liquid IgG fraction fraction Stabilizer: 0.2% BSA Preservative: 0.09% Sodium Azide Label: Fluorescein Isothiocyanate Isomer 1
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	FITC
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)