

Product datasheet for **SP2017HRP**

Porcine IgG F(c) specific Goat Polyclonal Antibody

Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgG F(c) specific Goat Polyclonal Antibody
Applications:	ELISA
Recommended Dilution:	ELISA: 1/10,000-1/100,000.
Reactivity:	Porcine
Host:	Goat
Immunogen:	Purified Porcine IgG.
Formulation:	PBS containing 0.09% ProClin™300 as preservative and 0.2% BSA as stabilizer Label: HRP State: Liquid purified IgG fraction Label: Horseradish Peroxidase
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	HRP
Storage:	Store undiluted at 2-8°C. DO NOT FREEZE!



[View online »](#)