

## Product datasheet for **SP2017P**

### Porcine IgG F(c) specific Goat Polyclonal Antibody

#### Product data:

Product Type:	Secondary Antibodies
Product Name:	Porcine IgG F(c) specific Goat Polyclonal Antibody
Applications:	ELISA, ID, WB
Recommended Dilution:	<b>ELISA:</b> 1/100-1/1000. <b>Western Blot:</b> 1/1000-1/30,000. <b>Immunodiffusion.</b>
Reactivity:	Porcine
Host:	Goat
Immunogen:	Purified Porcine IgG.
Formulation:	PBS containing 0.09% Sodium Azide as preservative State: Aff - Purified State: Liquid purified IgG fraction
Concentration:	lot specific
Purification:	Affinity Chromatography
Conjugation:	Unconjugated
Storage:	Store the antibody undiluted at 2-8°C for one month or (in aliquots) at -20°C for longer. Avoid repeated freezing and thawing.



[View online »](#)